


UNIQUE BRILLIANCE



ACUNEX[®] VARIOMAX


THE VARIFOCAL PREMIUM IOL
FOR BRILLIANT, NATURAL VISION

ACUNEX[®] VARIOMAX

EXCELLENT VISION DUE TO UNIQUE VARIFOCALITY


CONTINUOUS TRANSMISSION TECHNOLOGY

First IOL with varifocal effect and continuous surface profile for a homogeneous light distribution (high balanced effect)



SIGNIFICANT BETTER VISION

compared to a monofocal or trifocal IOL



WHAT DO EXPERTS SAY ABOUT THE ACUNEX® VARIOMAX:*



Prof. M. Rau (MD, Germany)

The best multifocal IOL for presbyopia correction I've ever used. Overall my patients are at most completely spectacle-free, even when reading the smallest fonts and are extremely satisfied. They have excellent visual results and almost no annoying photic effects.

Prof. J. Alio (MD, Spain)

The multifocal IOL that will never disappoint you! ACUNEX VarioMax is as easy to use as the simplest monofocal IOL, but offers the advantages of varifocal technology, with a continuous IOL power distribution. The excellent results at all distances lead to the best possible patient satisfaction.



* Source: Feedback of surgeons, www.teleon-surgical.com

WHAT DO PATIENTS SAY ABOUT THE ACUNEX® VARIOMAX:

Thomas W.,
age 51

After cataract surgery, my quality of life has increased enormously. Seeing so well exceeds all my expectations.

It feels like I am in my mid-20s again. Being able to see so clearly without glasses makes me incredibly happy.

Anna K.,
age 68

Ellen S.,
age 70


I do not need reading glasses anymore. I can even read the fine-print in my daily newspaper without any problems. Wonderful!

LIFESTYLE PATIENT SURVEY (N=595)**

How did your visual comfort in daily activities changed with the ACUNEX® VarioMax considering the following distances?


Far distance

[Driving at day / night, theater, cinema, hiking]




Intermediate distance

[Computer, gardening, cooking]



Near distance


[Reading, mobile phone, handicraft work]



■ Improved ■ Unchanged ■ Worse

** Archive data, based on study data with segmental MIOL

One-piece multifocal posterior chamber lens for optimised depth of focus with aspherical surface and blue light filter



Product	ACUNEX® VARIOMAX AN6VM	aspheric
Type	Foldable one-piece multifocal IOL for capsular bag fixation	multifocal
Optic Size	6.0 mm	blue light filter
Overall Length	12.5 mm	glistening-free
Haptic Angulation	0°	
Optic Design	Biconvex Sector-shaped nearvision segment - anterior: +3.0D Aspherical surface - posterior Spherical aberration neutral	
IOL Design	C-Loop haptic	
Material	Hydrophobic glistening-free acrylate UV absorbing with additional blue light filter	
Available Diopters	SE +10.0D to +30.0D (0.5D)	
Refractive Index	1.54 @ 35°	
A-Constant (nominal)	119.1	
Sterilisation	Gamma sterilisation	
Storage	Storage in 0.9 % saline solution	
Recommended Incision Size	≥ 2.0 mm	

BRILLIANT VISION WITHOUT COMPROMISE

ADVANTAGES OF THE ACUNEX® VARIOMAX

- **Natural vision** over all distances thanks to the unique optical design
- **Minimized photic phenomena** and **excellent contrast vision** day and night
- **High spectacle independence** for excellent vision in all everyday situations
- Premium care with **> 94% light transmission**: best in class of all multifocal IOLs
- **Maximum patient satisfaction**
- Positive word of mouth and thus **free advertising for your practice/clinic**

MANUFACTURER:

Teleon Surgical B.V.

✉ acunex@teleon-surgical.com

🌐 www.teleon-surgical.com